

# TRIOXINATOR™

## OPERATING INSTRUCTIONS

**READ ALL INSTRUCTIONS BEFORE USING**

**LIMITED MANUFACTURER WARRANTY  
- REGISTRATION REQUIRED -**



# READ ALL INSTRUCTIONS BEFORE USING

## DANGER

1. To avoid injury do not use the TRIOXINATOR™ if the glass probes are damaged.
2. Serious mechanical damage to the handle could damage the Li-ion battery. A damaged Li-ion battery could cause a fire.
3. A TRIOXINATOR™ with a mechanically damaged handle should be discarded immediately.
4. Do not use while bathing.
5. Do not place or store product where it can fall or be pulled into a tub or sink.
6. Do not place in or drop into water or other liquid.

**WARNING** - To reduce the risk of burns, electrocution, fire, or injury to persons.

1. Do not immerse the TRIOXINATOR™ in water or any liquid. Ensure that the TRIOXINATOR™ is kept away at all times from baths or basins filled with water. Also, do not store or place the TRIOXINATOR™ where it can fall into water or any liquid. If the TRIOXINATOR™ collects water or liquid it could be permanently damaged.

2. The TRIOXINATOR™ should never be left unattended while operating.
3. Close supervision is necessary when this product is used by, on, or near children or invalids.
4. The TRIOXINATOR™ is not for use by children.
5. Use this TRIOXINATOR™ only for its intended purpose as described in this manual. Never use a wall adaptor other than the supplied wall adaptor.
6. Never operate the TRIOXINATOR™ if it is not working properly or if it has been damaged externally by any means. Return the TRIOXINATOR™ to Follicles immediately for examination and repair.
7. Keep the cord of the wall adaptor away from heated surfaces.
8. Never use the TRIOXINATOR™ while sleeping or drowsy.
9. Never drop or insert any object or liquid into any opening of the TRIOXINATOR™.
10. Do not use the TRIOXINATOR™ where aerosol (spray) products are used or where oxygen is being administered.
11. Never try to open the TRIOXINATOR™.
12. Do not use the TRIOXINATOR™ directly on the eyes.
13. Consult your health practitioner if you are pregnant, or have a pace-maker, hearing device or any other medical device on your body, and wish to use the TRIOXINATOR™.

## **MAINTENANCE AND CLEANING**

1. Keep the glass probes clean for optimum operation. Regularly wipe the glass probes with a mild disinfectant. To avoid damage do not use excessive force when wiping the glass probes.
2. Make sure no cleaning liquid gets into the opening of the adaptor jack, as this could result in a malfunction.
3. Recharge the internal battery regularly for optimum service life.
4. If one or more of the glass probes are broken, send the TRIOXINATOR™ back to Follicles for inspection or repair. Mechanically damaged glass probes do not fall under warranty.
5. For optimum safety never try to repair the TRIOXINATOR™ by any means other than sending it back to Follicles for inspection or repair.

## **USER INSTRUCTIONS**

This product is intended for household use only. Its purpose is to enhance the appearance of hair and reduce the signs of hair loss.

### **NORMAL OPERATION**

1. The TRIOXINATOR™ only operates when the wall adaptor is disconnected.
2. There are 5 output settings that can be selected using the push button.
3. Activating the push button while on position 5 terminates operation.

4. The letter “E” indicates an empty battery.
5. The TRIOXINATOR™ stops operating after 15 Minutes.

### **BATTERY CHARGING**

1. Connect the wall adaptor (110-220V 50/60Hz) to the TRIOXINATOR™.
2. The letter “C” indicates the battery is charging.
3. The letter “F” indicates a fully charged battery.
4. If the TRIOXINATOR™ remains uncharged for more than one year it might go into the pre-charge mode initially indicated by a flashing “E” on the display.
5. The TRIOXINATOR™ will not charge at temperatures below 0°C and above 45°C.
6. An empty battery will become fully charged in less than 4 hrs.
7. A new fully charged battery lasts for at least 8 sessions of 15 min. each.
8. A new TRIOXINATOR™ needs charging before use.

### **OPERATING CONSIDERATIONS**

1. The suggested massaging time for the affected area is 15 minutes every second day. More frequent or longer massaging does not necessarily increase effectiveness.
2. Adjust the power setting for a comfortable tingling effect. Going beyond this point may or may not enhance results.

3. Proper operation requires that the palm of the hand firmly touch the top of the handle. The top of the handle is on the same side as the push button.
4. When massaging the scalp make sure the glass probes touch the skin as much as possible.
5. When massaging areas with existing hair make sure that the glass probes are touching the skin as much as possible.
6. Gently and evenly move the TRIOXINATOR™ over the entire area of the scalp to ensure that the scalp area is massaged evenly for 15 minutes. Do not apply extra pressure on the skin. Regular cleaning of the glass probes facilitates smooth movement over the skin.
7. Do not use any product on the hair when massaging the scalp with the TRIOXINATOR™.
8. The hair must be clean and dry when using the TRIOXINATOR™.

### **ADDITIONAL INFORMATION**

In the highly unlikely situation that dirt gets into the adaptor jack resulting in malfunctioning of the unit, the following procedure can be helpful:

1. Turn off TRIOXINATOR™ – ensure battery is not empty
2. Plug in live wall adaptor.
3. Press and hold the push button until the vertical line appears on display.
4. Release push button while vertical line is visible.

### **RETURNING TO ORIGINAL CONDITION**

1. Turn off TRIOXINATOR™ – ensure battery is not empty.
2. Plug in Wall adaptor.
3. Press and hold the push button until the vertical line appears on display.
4. Release push button after vertical line disappears.

## **SAVE THESE INSTRUCTIONS**

**See next page for  
Limited Manufacturer  
Warranty Information**

## IMPORTANT - PLEASE READ! LIMITED MANUFACTURER WARRANTY

This warranty extends only to the original purchaser of the TRIOXINATOR™ and is limited to the repair or replacement of the device. Online registration of warranty at [www.follicleshair.com](http://www.follicleshair.com) is required in order to activate warranty.

Go to <http://www.follicleshair.com> and click on the “Warranty” Link to activate your warranty. You will be required to provide the serial number from your unit. The Serial Number (S/N) is contained on this label.



**NOTE: Proof of purchase in the form of a bill of sale or receipted invoice which is evidence that the unit is within the warranty period must be presented to obtain warranty service.**

## IMPORTANT - PLEASE READ! LIMITED MANUFACTURER WARRANTY

Extended warranty available on-line at the time of registration. Terms of which are provided on-line only at <http://www.follicleshair.com>

FOLLICLES WORLDWIDE MARKETING GROUP INC. (“FOLLICLES”) Limited Manufacturer Warranty warrants this Product (including any accessories) against defects in material or workmanship as follows:

**LABOR:** For a period of one (1) year for the original purchaser from the date of purchase, if this Product is determined to be defective, FOLLICLES will repair or replace the Product, at its option, at no charge, or pay the labor charges to any FOLLICLES authorized service facility. After the Warranty Period, you must pay for all labor charges.

**PARTS:** FOLLICLES will supply, at no charge, new or rebuilt replacements in exchange for defective parts for a period of one (1) year from the date of purchase.

**SHIPPING:** This warranty does not include shipping fees.

### EXCLUSIONS:

The following is not covered by the Limited Manufacturer Warranty:

**NORMAL WEAR and TEAR:** Periodic maintenance, repair and replacement of parts due to normal wear and tear are excluded from coverage.

**BATTERIES:** Only batteries whose fully charged capacity falls below 80% of their rated capacity within the warranty duration and batteries that leak are covered by this limited warranty.

**ABUSE and MISUSE:** Defects or damage that result from: (a) broken or cracked glass tubes, improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; (b) contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; (c) use of the Products or Accessories for commercial purposes or subjecting the Product or Accessory to abnormal usage or conditions; or (d) damage due to connection to improper voltage supply, or attempted repair by anyone other than a facility authorized by FOLLICLES to service the Product; or (e) other acts which are not the fault of FOLLICLES, including acts of God, are excluded from coverage.

**USE of NON-FOLLICLES PRODUCTS & ACCESSORIES:** Defects or damage that result from the use of Non-FOLLICLES branded or certified Products, Accessories or other peripheral equipment are excluded from coverage.

**UNAUTHORIZED SERVICE or MODIFICATION:** Defects or damages resulting from service, testing, adjustment, maintenance, alteration, including without limitation, modification in any way by someone other than FOLLICLES, or its authorized service centers, are excluded from coverage.

**ALTERED PRODUCTS:** Products or Accessories with (a) serial numbers or date tags that have been removed, altered or obliterated; (b) broken or tampered seals; (c) mismatched board serial numbers; or (d) nonconforming or non-FOLLICLES housings or parts, are excluded from coverage.

**EFFICACY:** This limited manufacturer warranty does not cover customer instruction or product efficacy.

**REPAIR OR REPLACEMENT AS PROVIDED UNDER THIS WARRANTY IS THE EXCLUSIVE REMEDY OF THE CONSUMER. FOLLICLES SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES FOR BREACH OF ANY EXPRESS OR IMPLIED WARRANTY ON THIS PRODUCT.**

Follicles Worldwide Marketing Group Inc.  
260-4255 Arbutus Street  
Vancouver BC V6J 4R1  
Canada

Telephone: 1-866-611-5678  
Vancouver BC - 604-732-1183

Email: [info@follicleshair.com](mailto:info@follicleshair.com)

<http://www.follicleshair.com>



Made in Canada

